



TENEX
HEALTH

RESTORE FUNCTION WITH
Tenex Health TX®

Treat diseased soft tissue and restore musculoskeletal
function without invasive surgery

Join us at Booth #0000

Conference Name 2019 | City

Convention Center Name

Exhibit Dates x-xx

Visit our booth and demo of a new level of cutting with the *TX2 MicroTip* for Shoulder and Hip Tendinosis

You are cordially invited to a hands-on demonstration of the *Tenex Health TX System®* and *TX MicroTip*. Over 85,000 procedures performed since 2012.


During the demonstration you will have a hands-on opportunity to use the innovative *Tenex Health TX* technology to:

- Identify and simultaneously perform specific cutting and removal of pain-generating pathologic soft tissue
- Experience how the *TX2 MicroTip* with longer needle length (2 inches) can minimally invasively cut and remove previously unreachable tendon tissue, such as in the shoulder and hip
- Visualize how the *TX1* and *TX2 MicroTip* definitively treats diseased tendon tissue without disturbing surrounding healthy tendon, while stimulating healthy tissue growth



Common applications for *Tenex Health TX®*

- Shoulder Tendinosis
- Lateral or Medial Epicondylitis
- Gluteal Tendinosis of the Hip
- Patellar Tendinosis
- Achilles Tendinosis
- Plantar Fasciosis



Calcific Deposit Tendon Tissue



TX2 MicroTip (2 inches)

Tenex Health TX® for Percutaneous Tenotomy and Fasciotomy – is based on technology developed in collaboration with the Mayo Clinic. Tenex Health TX combines conventional ultrasound imaging for visualization with the advanced TX MicroTip for rapid and precise cutting and removal of diseased soft tissue.

Tenex Health, Inc.
26902 Vista Terrace, Lake Forest, CA 92630
Corporate: 949.454.7500
www.TenexHealth.com
(MKT-288, Rev. A)